

Certificate of Analysis

Catalogue Number

KL-02-00038

Lot Number

KLAH00117

Product Name

Adapalene Related Compound A

Synonyms

Methyl 6-bromonaphthalene-2-carboxylate

CAS Number

33626-98-1

Molecular Formula
 $C_{12}H_9BrO_2$
Molecular Weight


265.1 g/mol

Mfg. Date

06-10-2025

Re-Test Date

04-10-2028

Structure

Shipping Condition

All Product are stable to be shipped at room temperature, unless otherwise specified.


Long Term Storage Conditions

2°C ~ 8°C Protected from light and Sealed

Analytical Information

Test	Specifications	Results
Appearance	White to off white solid	Off-white crystalline powder
Solubility	Soluble in Chloroform	Conforms
Purity by HPLC	Not less than 90.00%	99.89%
Identification by 1H NMR	Conforms to structure	Complies
Identification by Mass	Conforms Molecular Weight	Complies
Identification by IR	Conform functional groups	Complies
Weight Content By TGA	Not specified	0.031%
Residue of Ignition By TGA	Not specified	0.000%
Potency	Not specified	99.86%

Prepared by



Reviewed & Approved by



Note: Use for research & analytical purpose only. Not for human consumption

Kaaris Research Labs Private Limited, Shed No. 2B at Plot No.22, Road No. 17, IDA, Nacharam, Hyderabad-500 076

sales@kaarislabs.com, www.kaarislabs.com An ISO 9001-2015 Certified Company