

Certificate of Analysis

Catalogue Number

KL-05-00035

Lot Number

KLAJ00206

Product Name

Bromobenzene-d5

Synonyms

1-bromobenzene-2,3,4,5,6-d5

CAS Number

4165-57-5

Molecular Formula
 C_6D_5Br
Molecular Weight


162.04 g/mol

Mfg. Date

30-07-2025

Re-Test Date

25-01-2026

Structure

Shipping Condition

All Product are stable to be shipped at room temperature, unless otherwise specified.


Long Term Storage Conditions

2°C ~ 8°C Protected from light and Sealed

Analytical Information

Test	Specifications	Results
Appearance	Clear Colourless Oil	Clear Colourless Oil
Purity by GC (% area)	Not less than 99.00%	99.90%
Isotopic Purity by QNMR	Not less than 99.00%	99.10%
Water content by KF (% w/v)	Not more than 0.5	0.03
Identification by 1H NMR	Conforms to structure	Complies

Prepared by



Reviewed & Approved by



Note: Use for research & analytical purpose only. Not for human consumption

Kaaris Research Labs Private Limited, Shed No. 2B at Plot No.22, Road No. 17, IDA, Nacharam, Hyderabad-500 076

sales@kaarislabs.com, www.kaarislabs.com An ISO 9001-2015 Certified Company