

Certificate of Analysis

Catalogue Number

KL-05-00033

Lot Number

KLAK00282

Structure
Product Name

Iodomethane-D3

Synonyms

iodomethane-d3

CAS Number

865-50-9

Molecular Formula

 CD₃I

Molecular Weight


144.96 g/mol

Mfg. Date

01-03-2025

Re-Test Date

01-03-2027


Shipping Condition

All Product are stable to be shipped at room temperature, unless otherwise specified.


Long Term Storage Conditions

2°C ~ 8°C Protected from light and Sealed

Analytical Information

Test	Specifications	Results
Appearance	Colorless or brown or pink Liquid	Colorless or brown or pink Liquid
Purity by GC	Not less than 99.00%	99.90%
Isotopic Purity	Not less than 99.50%	99.50%
Identification by ¹ HNMR	Conforms to structure	Complies
Moisture Content (Karl Fisher)	Max. 300 ppm	200 ppm

Prepared by



Reviewed & Approved by



Note: Use for research & analytical purpose only. Not for human consumption

Kaaris Research Labs Private Limited, Shed No. 2B at Plot No.22, Road No. 17, IDA, Nacharam, Hyderabad-500 076

sales@kaarislabs.com, www.kaarislabs.com An ISO 9001-2015 Certified Company