

Certificate of Analysis

Catalogue Number

KL-06-00311

Lot Number

KLP0050AJ01

Product Name

N-Nitroso Risperidone Difluorobenzoyl Impurity

Synonyms

(2,4-difluorophenyl)(1-nitrosopiperidin-4-yl)methanone

CAS Number

NA

Molecular Formula
 $C_{12}H_{12}F_2N_2O_2$
Molecular Weight


254.24 g/mol

Mfg. Date

12-10-2025

Re-Test Date

10-10-2028

Structure

Shipping Condition

All Product are stable to be shipped at room temperature, unless otherwise specified.


Long Term Storage Conditions

2°C ~ 8°C Protected from light and Sealed

Analytical Information

Test	Specifications	Results
Appearance	Pale Yellow Solid	Pale Yellow Solid
Solubility	Soluble in DMSO, Methanol	Conforms
Purity by HPLC	Not less than 90.00%	97.69%
Identification by ^1H NMR	Conforms to structure	Complies
Identification by Mass	Conforms Molecular Weight	Complies
Identification by IR	Conform functional groups	Complies
Weight Content By TGA	Not specified	0.318%
Residue of Ignition By TGA	Not specified	0.000%
Potency	Not specified	97.38%

Prepared by



Reviewed & Approved by



Note: Use for research & analytical purpose only. Not for human consumption

Kaaris Research Labs Private Limited, Shed No. 2B at Plot No.22, Road No. 17, IDA, Nacharam, Hyderabad-500 076

sales@kaarislabs.com, www.kaarislabs.com An ISO 9001-2015 Certified Company