

## Certificate of Analysis

**Catalogue Number**

KL-05-00046

**Lot Number**

KLP0076AJ01

**Product Name**

Phenyl-d5-boronic acid

**Synonyms**

(phenyl-d5)boronic acid

**CAS Number**

215527-70-1

**Molecular Formula**
 $C_6H_5D_5BO_2$ 
**Molecular Weight**


126.96 g/mol

**Mfg. Date**

08-11-2025

**Re-Test Date**

07-11-2027

**Structure**

**Shipping Condition**

All Product are stable to be shipped at room temperature, unless otherwise specified.


**Long Term Storage Conditions**

2°C ~ 8°C Protected from light and Sealed

**Analytical Information**

| Test                           | Specifications           | Results                                             |
|--------------------------------|--------------------------|-----------------------------------------------------|
| Appearance                     | White to Off-White Solid | White to Off-White Solid                            |
| Purity by HPLC                 | Not less than 95.00%     | 99.00%                                              |
| Isotopic enrichment (by Q-NMR) | Not less than 95%        | 98.95% (Internal standards is 1,4-dimethoxybenzene) |
| Identification by $^1H$ NMR    | Conforms to structure    | Complies                                            |

Prepared by



Reviewed &amp; Approved by



Note: Use for research &amp; analytical purpose only. Not for human consumption

Kaaris Research Labs Private Limited, Shed No. 2B at Plot No.22, Road No. 17, IDA, Nacharam, Hyderabad-500 076

sales@kaarislabs.com, www.kaarislabs.com An ISO 9001-2015 Certified Company